

BASAL MEDIUM 2

	<u>1 litre</u>	<u>20 ml</u>
Phosphate ammonia base (BM2 10X)	100 ml	2 ml
H ₂ O	900	18
Carbon source: 40% glucose	10	0.2
50 mM MgSO ₄	10	0.2
10 mM FeSO ₄	1	0.02

Wash all glassware at least three times with dH₂O.

Make and autoclave all ingredients separately and add aseptically.

40% glucose becomes 0.4% glucose.

Note: When making 10mM FeSO₄, add a few mls of concentrated HCl and the iron will stay in solution. (Red color will almost disappear – if not be sure to shake well before adding to media)

BM2 10X

	<u>1 litre</u>
0.07 M (NH ₄) ₂ SO ₄	9.25 g
0.4 M K ₂ HPO ₄	69.7
0.22 KH ₂ PO ₄	29.9
pH 7.0	

***ALTERNATE CARBON SOURCES*:**

	<u>1 litre</u>
2 M K succinate, pH 7.0	10 ml
2 M K gluconate, pH 7.0	10
2 M K pyruvate, pH 8.0	10

Make each from their respective acids (succinic acid, gluconic acid, pyruvic acid) and pH using KOH pellets. You will need a fair amount of KOH to reach a neutral pH.

2 M solutions become 0.02M.

You can autoclave glucose, K gluconate, and K succinate but filter sterilize K pyruvate. Check pH of older stocks.