

M9 MINIMAL MEDIUM

REFERENCE: Molecular Cloning, A Laboratory Manual 1st ed. Maniatis. 1982. p.68.

STOCK SOLUTIONS:

- | | | | |
|----|--------------|----------------------------------|--|
| 1. | 10X M9 SALTS | Na ₂ HPO ₄ | 60 g |
| | | KH ₂ PO ₄ | 30 g |
| | | NaCl | 5 g |
| | | NH ₄ Cl | 10 g |
| | | dH ₂ O | to 1 litre (pH before adding
all the water) |

Adjust to pH 7.4 with 10 M NaOH or NaOH pellets.
Autoclave to sterilize.

- | | | |
|----|-----------------------|-------------------|
| 2. | 1 M MgSO ₄ | Autoclave. |
| 3. | 1 M CaCl ₂ | Autoclave. |
| 4. | 40% (w/v) glucose | Filter sterilize. |

FOR BROTH:

Autoclave 900 mls dH₂O.
Let cool then aseptically add:

100 mls M9 salts
2 mls 1 M MgSO₄
0.1 ml 1 M CaCl₂
5 mls 40% glucose

Note: Do Not add CaCl₂ for Pseudomonas

FOR AGAR:

Autoclave 900 mls dH₂O
20 g agar

Let cool to 50°C then aseptically add ingredients as for broth.

- NOTES:**
- CaCl₂ will precipitate if added when media is too hot.
 - You can autoclave dH₂O and 10X salts together when making the broth, then add the rest when cool.