

DEFINED PHOSPHATE MINIMAL MEDIA (KEITH'S)

For *P. aeruginosa* and other bacteria

REFERENCE: J. Bacteriology. 150:730, 1982.

		1 litre
1 M Hepes (pH 7.4, use HCl or 5 N NaOH)		100 ml
10X (NH ₄) ₂ SO ₄ (=70 mM; 10X = 9.25 g/l)		100
2 M K succinate (pH to 7.2 with solid KOH)		10
0.1 M MgSO ₄		5
10 mM FeSO ₄		1
Ion solution		1
BM2 10X (no (NH ₄) ₂ SO ₄)	Pi suf. = 1.5 mls	Pi def. = 0.3 mls
dH ₂ O	Pi suf. = 781.5 mls	Pi def. = 783 mls

Autoclave dH₂O, Hepes + (NH₄)₂SO₄ together, others separately (except ion solution which must be chloroform sterilized). Add others after autoclaving.

ION SOLUTION

REFERENCE: Hancock, Raffle and Nicas. AAC. 19(5):777 – 785, 1981.

1.78 mM	FeCl ₃ -5H ₂ O	48 mg
1.62 mM	MnCl ₂ -4H ₂ O	32 mg
2.45 mM	CaCl ₂ -2H ₂ O	36 mg
13.91 mM	ZnSO ₄ -7H ₂ O	400 mg
4.69 mM	H ₃ BO ₃	29 mg
0.67 mM	CsCl-6H ₂ O	16 mg

In 100 mls dH₂O

Sterilize by adding 5 ml of CHCl₃ and shaking.